



Instructions For Use

VirSNIp SARS-CoV-2 Spike del H69/V70

530

Cat.-No. 53-0781-96

Roche SAP n° 09 405 674 001

Kit with reagents for 96 PCR reactions 20 µl for genotyping of SARS-CoV-2 RNA [lyophilized]

1. Content, Storage and Expiry

- 1 Vial yellow cap 96 reactions CoV (lyophilized)
1 Mixed positive control wt/69,70del

Storage at Arrival:

Store cooled or at ambient temperature
Do not freeze the lyophilized reagents.

- Kits are stable for one year after production (store 4°C to 25°C in the dark). See lot-specific expiry date.
• Reconstituted reagents are stable for two weeks if stored protected from light and cooled (2°C to 8°C).
• Dissolved reagent can be stored long-term if frozen (-15°C to -25°C). Avoid multiple freeze-thaw cycles.
• Reconstituted positive controls must be stored frozen. Minimize multiple freeze-thaw cycles.

2. Additional Reagents required

LightCycler® Multiplex RNA Virus Master
or 1-step RT polymerase

Roche Cat.-No. 06 754 155 001
90-9999-96

3. Introduction

The SARS-CoV-2 genome was published 11.1.20 (Genbank MN908947). Hundreds of thousand isolates have been sequenced. The deletion 69/70 was first reported für Alpha. Omicron bears an additional 67V.

Table with 10 columns: Spike Prot. Variation, Genetic Variation, UK B.1.1.7, UK B.1.525, ZA B.1.351, Brazil P.1, Brazil P.2, Omicron B.1.1.529, Function, Effect, Single assay. Rows include variations like A67V, del HV69/70, K417N, N439K, Y453F, E484K, N501Y, A570D, D614G, P681H, V1176F.

Proof of the deletion alone is no evidence for the UK strain. Variant is frequent (eg. Berchtesgarden)

4. Description

A 119 (del: 113) bp long fragment is amplified and analyzed with a melting curve using a deletion-specific detection probe. With the Roche polymerase the amplification of the wild type variant is barely visible.

5. Specification

Sensitivity better than 50 copies viral RNA UK strain B.1.1.7 (50% signal compared to 4,000 copies).

6. Sample Material and Extraction

Coronaviruses affect normally the lower respiratory system, but SARS-CoV-2 is found also in the upper part. Typical clinical samples are throat or nasopharyngeal swabs, sputum, saliva, or gargle solution. Product tested with heat-treated gargle solution. For RNA extraction see manufacturer's kit instructions.

7. Material Safety Data (MSDS)

This product is not hazardous (according to regulation (EC) No 1272/2008), not toxic, not IATA-restricted. Not from human, animal or plant origin. Product contains synthetic oligonucleotide primers and probes. According to OSHA 29CFR1910.1200, Australia [NOHSC:1005, 1008 (1999)] and EU Directives (EC) No 1907/2006 and (EC) No 2015/830 any products which do not contain more than 1% of a component classified as hazardous or classified as carcinogenic do not require a MSDS.



8. Instructions for Use

Instruction for Roche 480 instruments. Capillary LightCycler[®], LightCycler[®] 96, MyGo and BioRad CFX96 instruments give similar results. For other instruments use the SYBR Green melting option.

8.1. Programming Roche 480 Instruments

Detection Format 530 Channel **Set Quant Factor 10, Max Integration Time 1 sec**

LightCycler[®] 480 Instrument: 483-533

LightCycler[®] 480 II Instrument: 465-510

cobas z 480 Analyzer (open channel): 465-510

Program Step:	RT Step	Denaturation	Cycling			Cooling
Parameter						
Analysis Mode	None	None	Quantification mode			None
Cycles	1	1	40-45			1
Target [°C]	55	95	95	60	72*	40
Hold [hh:mm:ss]	00:05:00	00:05:00	00:00:05	00:00:15	00:00:15	00:00:30
Ramp Rate [°C/s] 96	4.4	4.4	4.4	2.2	4.4	1.5
Ramp Rate [°C/s] 384	4.6	4.6	4.6	2.4	4.6	2.0
Acquisition Mode	None	None	None	Single	None	None

* 72°C step can be skipped. 95°C can be cut to 3 s, 60°C to 12 s. RT and Den to 3 min (total time 45 min) Table 1

8.1.1. Melting Analysis (may be added or programmed as second run)

Detection Format Hydrolysis Probe or SimpleProbe

LightCycler[®] 480 Instrument: 483-533

LightCycler[®] 480 II Instrument: 465-510

cobas z 480 Analyzer (open channel): 465-510

Program Step:	Melting			Cooling
Parameter				
Analysis Mode	Melting Curves mode			None
Cycles	1			
Target [°C]	95	40	75	40
Hold [hh:mm:ss]	00:00:30	00:02:00	00:00:00	00:00:30
Ramp Rate [°C/s]	4.4	1,5	-	1.5
Acquisition Mode	-	-	Continuous	
Acquisitions [per °C]	-	-	3**	None

Table 2

**Melting slope shall be 0.19 to 0.29°C per second. If reading more channels reduce the number of acquisitions/sec.

8.2. Experimental Protocol

- **Sample material:** Use aqueous nucleic acid preparations
- **Negative control:** Always run at least one no-template control (NTC) - replace the template NA with water.
- **Positive control:** Run a positive control - replace the template NA with the provided Positive Control.

For an increased sensitivity use 10 µl nucleic acid per 20 µl reaction.

Product tested for 10 µl reaction volume (192 reactions).

8.2.1. Preparation of Parameter-Specific Reagents (PSR, 96 reactions):

The reagent vial with the **yellow** cap contains the primers and probe to run 96+ PCR reactions.

Check for the orange pellet, then **add 50 µl** PCR-grade water, mix (vortex) and spin down. For robotic pipetting the volume can be extended to 55 µl (signals will decrease by 10-20%).

► **Use 0.5 µl** reagent per 20 µl PCR reaction.

8.2.2. Preparation of the Positive Control

Add 160 µl PCR-grade water to the vial with the **black** cap, if using 10 µl sample volume add 320 µl. Mix by pipetting up and down 10 times. If vortexing spin down to collect the solution. Store frozen.

Notes: Opening this vial may cause contamination of the workspace. Pulse spin vial prior to opening.

► Use 5 µl positive control for a 20 µl PCR reaction (10 µl if using 10 µl sample volume).

8.2.3. Preparation of the Reaction Mix

Multiply volumes by the number of reactions plus controls and one reserve, prepare in a cooled tube:

LightCycler® Multiplex RNA Virus Master		1-step RT Polymerase 90-9999-96	
10.4 µl	Water, PCR-grade	Water, PCR-grade	4.5 µl
0.5 µl	This reagent (PSR)	This reagent (PSR)	0.5 µl
4.0 µl	Roche Mastermix	TIB Mastermix	10.0 µl
0.1 µl	RT Enzym	-	-
15.0 µl		15.0 µl	

Table 3

Mix gently, spin down and transfer 15 µl (10 µl) per well.

Add 5 µl (10 µl) of sample or control to each well for a final reaction volume of 20 µl. Seal plate and centrifuge.

Start run

9. Typical Results obtained with Roche Polymerase




Figure 1. Blue SARS RNA, Green wt, Red deletion **Left** Deletion is clearly visible in the amplification. **Center:** Melting curves. **Right** Wt (and Delta) has a melting temperature of 56.5 (± 2)°C, del69/70 (Alpha) has a Tm of 64.0 (± 2)°C, 67V del69/70 (Omicron) has a Tm of 58.0 (± 2)°C

10. Reading the Results

Use the Second Derivative Maximum method (Automated (F'' max). View results in the 530 channel. The negative control (NTC) must show no signal. For the melting curve analysis use 'Tm calling'.

Channel 530 Amplification	Channel 530 Melting analysis	Channel 530 NTC Control	Result
Not relevant	Not relevant	Negative / no peak	No virus amplified / not detectable
Invisible	Tm ~ 56°C*	Negative	for example C21771T (no AA change)
Invisible or low sig.	Tm ~ 58°C*	Negative	SARS Spike wild type (not UK variant)
Visible	Tm ~ 64°C*	Negative	SARS Spike del69/70 (UK or cluster V)
Not relevant	Not relevant	Positive	Contamination Repeat experiment

Tm values shift depending on the instrument, speed of heating, mastermix, salt contents and detection format.

* Temperatures with 1step RT pol. 90-9999-96 are 3-4°C higher.

Single peak with lower Tm values are an indication for the presence of another mutation in the probe region.

The assay detects the genetic situation and not a strain; the correlation to a reference strain describes the most likely assignment for isolates reported as VOI and VOC.

11. References

Genomic characterisation of emergent SARS-CoV-2 lineage in UK defined by novel set of spike mutations. Rambaut et al., 2020 www.ecdc.europa.eu/en/publications-data/threat-assessment-brief-rapid-increase-sars-cov-2-variant-united-kingdom

Circulating SARS-CoV-2 spike var. N439K maintains fitness while evading antibody-mediated immunity. Thomson et al., 2020

Mutations in SARS-CoV-2 spike protein and RNA pol. are associated with COVID-19 mortality risk. Hahn et al., 2020 www.rki.de/DE/Content/InfAZ/N/Neuartiges_Coronavirus/DESH/Bericht_VOC_05022021.pdf?__blob=publicationFile

12. Multiplex PCR Compatibility


This SNP assay can be combined with 51-0776-96 SARS E+N and either EAV spiked extraction control or UBC human mRNA extraction control or with the complete kit 60-0770-96 Sarbecovirus E+N+UBC.

Multiplex PCR and Instrument Compatibility						480 II	z 480	LC96	LC2.0	Nano
Color Comp 40-0320 mandatory only for Multiplex PCR using more channels										
500	530	580	610	640	660					
	SNP					X	X	X	X	X
SarbecoV	SNP	UBC mRNA				X	X	X		
SarbecoV	SNP				UBC	X	X	X		
SarbecoV	SNP				EAV	X	X	X		
	SNP	SARS N	SARS E		UBC	X	X	X		

Table 3

13. Version History

V201227	Release version	2020-12-22
V210122	10. Tm values for 1-step RT polymerase 90-9999-96	2021-01-22
V210123	8.2.3 Instructions for 1-step RT polymerase 90-9999-96	2021-02-04
V210404	Variant table 9. Tm for different variants	2021-04-20
V211202	Omicron variant has a Tm 2°C higher than wild type or Delta	2021-12-02

Certificate of Analysis (CoA)						
Lot n° 4996 Expiry : YYYY-MM-DD						
	wt	del69/70	67V d69/70		PC	passed
Tm range Measured	56-59°C	63-66°C	58-59°C	Cp range	22-30	✓
Signal level Measured	2-10	2-10	2-10			✓
Negatives	10/10					✓
<p>Note: Cp (crossing point) values collected with pDNA (single target PCR). Fluorescence (FL) levels depend on instrument settings and may vary. The Cp values will vary from instrument to instrument by up to 2 cycles, while the distance between two dilution steps should be relative constant (ΔCp).</p>						
DOM (manufactured): YYYY-MM-DD			QC Acceptance: YYYY-MM-DD			
We, the undersigned, certify that the product designated above has been obtained in accordance with the rules of production and quality control.						
Name(s) :						
Name1			Name2			

TIB MOLBIOL Syntheselabor GmbH | Eresburgstr. 22-23 | D-12103 Berlin | Germany
 Tel. +49 30 78 79 94 55 | FAX +49 78 79 94 99 | dna@tib-molbiol.de | WWW.TIB-MOLBIOL.COM
 HS code 3822 0000 | EORI DE 4806433 | Registry Court Berlin Charlottenburg HRB 93163 B